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Child/Adult mask - PVC

### COMPONENTS-(FIG. A)

- |                             |                         |
|-----------------------------|-------------------------|
| 1. Compressor housing       | 10. Adult mask          |
| 2. Air outlet               | 11. Pediatric mask      |
| 3. Battery level indicators | 12. Replacement filters |
| 4. Filter compartment       | 13. Nosepiece           |
| 5. Micro USB port           | 14. Micro USB cable     |
| 6. ON/OFF button            | 15. Power adapter       |
| 7. Nebulizer                | 16. Storage bag         |
| 8. Mouthpiece               |                         |
| 9. Air tube                 |                         |

**IMPORTANT:**

Read the instructions carefully before use.

**INTENDED USE**

AirForce Mini is an aerosol therapy system indicated for home use.

**WARNINGS**

1. Use the device only as described in this manual, therefore as an aerosol therapy system, and follow the directions provided by your physician. Any use different from the intended one is to be considered improper and hence dangerous; the manufacturer cannot be held liable for any damage caused by improper, incorrect and/or unreasonable use, or if the equipment is connected to electrical installations which do not comply with current safety regulations.
2. Keep this manual for future reference.
3. Do not use the device in the presence of any anaesthetic mixture inflammable with oxygen or nitrogen protoxide.
4. The correct operation of the device may be affected by electromagnetic interferences that exceed the limits indicated by the European standards in force. In case this device interferes with other electrical devices, move it and plug it to a different power socket.
5. In case of failure and/or malfunction of the device, read the "TROUBLESHOOTING" section. Do not tamper with or open the compressor's housing.
6. For repairs address only to a technical service centre authorised by the manufacturer and always request the use of

original spare parts. Failure to comply with the above mentioned indications can compromise the safety of the device.

7. Comply with the safety regulations concerning the electrical devices and in particular:
  - use only original accessories and components;
  - never immerse the device in water;
  - never wet the device, it is not protected against water penetration;
  - never touch the device with wet or moist hands;
  - do not leave the device exposed to atmospheric agents;
  - the use of this device by children and disabled people always requires the close supervision by an adult with full mental faculties;
  - do not pull the power cord or the device itself to unplug it from the power socket;
8. The device operates on batteries. You can also use the device plugged to the mains.
9. Before connecting the power plug, ensure that the electrical data shown on the rating plate on the bottom of the unit, matches the nominal data of the mains.
10. If the plug supplied with the appliance is not compatible with the mains electricity socket, contact qualified staff to replace the plug with a suitable one. Mainly, it is not advisable to use adaptors, either single or multiple, and/or extension leads. If the use thereof is necessary, you must use types conform to safety standards, being however careful not to exceed the maximum specified power limits that are written on the adaptors and the extension leads.

11. Never leave the appliance plugged when not used: unplug the device from the mains outlet when the device is not in use or has fully charged.
12. The installation has to be carried out by qualified personnel following the instructions of the manufacturer. Incorrect installation can cause damage to people, animals or property, for which the manufacturer cannot be deemed liable.
13. The USB cable and the power supply adapter of this device must not be replaced by the user. In case of damage, have them replaced by a service centre authorised by the manufacturer.
14. It is recommended to fully extend the power cord, throughout its entire length, to avoid dangerous overheating.
15. Before performing any cleaning and/or maintenance operation, switch off the device and unplug it.
16. Some parts of the device are so small that they can be swallowed by children; keep the device out of the reach of children.
17. Should you decide not to use the device, it is recommended to dispose it of in accordance with the current regulations.
18. Make sure to:
  - use this device only with medicines prescribed by your doctor;
  - carry out the treatment only using the accessory indicated by your doctor depending on the pathology;
  - use the nosepiece accessory only if expressly indicated by your doctor and be careful to NEVER introduce the

bifurcations in the nose, but only bring them as close as possible to it.

19. Check in the medicine package leaflet for any contraindications for use with common aerosol therapy systems.
20. Do not place the device so that it is difficult to unplug.
21. To prevent strangulation and suffocation, keep the power cord and the air hose out of the reach of children.
- 22. The accessories are to be used for a single patient. It is not recommended to use them on several patients.**

## HOW TO USE THE DEVICE

The device must be checked before each use in order to identify any malfunctions and/or damage caused during transportation and/or storage. During inhalation, sit upright and relaxed to avoid compressing the airways and impairing treatment effectiveness.

1. After unpacking the device, check it for visible damages or defects; pay particular attention to cracks in the plastic housing, which may expose electrical components. Check if the accessories are intact.
2. Before using the device, proceed with the cleaning operations as described in the chapter «CLEANING AND MAINTENANCE».

3. Open the nebulizer by turning counterclockwise the top (Pic. B).

4. Make sure that the medicine conduction cone is properly fitted on the air conduction cone inside the nebulizer (Pic. C).

5. Pour into the nebuliser the amount of medicine indicated by your doctor (Pic. D).

6. Close the nebulizer by turning clockwise the two parts making sure that they are properly sealed (Pic. E)

7. Connect one end of the air tube to the dedicated outlet at the bottom of the nebulizer (Pic. F) and the other end to the air outlet (Pic. G) on the device

8. If you need to use the mask, insert it directly onto the nebulizer (Pic. H).

9. To start treatment, press the I/O switch and make sure the red light around the power button is steadily lit (PIC.I)

10. The battery inside the device is indicated by the blue LED below the battery icon (PIC.L)

11. If the battery is empty, the blue LEDs are not lit and the orange light around the I/O switch flashes (PIC.M)

12. When the battery is empty, plug the device to the mains outlet and wait a few minutes before you switch it on (PIC.N)

13. While charging, the blue LED flashes (PIC.O).

14. In case of prolonged non-use, we recommend that you charge the battery completely every 3 months to ensure maximum performance of the same.

 **CAUTION! The device is designed for intermittent use (30 min. on - 30 min. off)**

15. Inhale the aerosol solution using the prescribed accessory.
16. Once you have finished the treatment, switch the appliance off by pressing the I/O switch.
17. Clean the nebulizer and its accessories following the steps in the «CLEANING AND MAINTENANCE» section.
18. The device requires no calibration.
19. Tampering with the device is absolutely forbidden.

## **CLEANING AND MAINTENANCE**

The device must be cleaned with a soft dry cloth and non-abrasive cleaners.

 **CAUTION!**

During cleaning operations, make sure the internal parts of the device are not in contact with liquids and that the power plug is disconnected.

### **Cleaning and disinfection of the accessories**

 Follow the directions on cleaning and disinfection of the accessories as they are critical to product performance and therapy success.

#### **1. Upon first use and at the end of each application:**

Disassemble the nebulizer by turning the top counterclockwise and remove the medication conduction cone. Wash the components of the disassembled nebulizer, the mouthpiece and nosepiece with tap water. Finally, dip them in boiling water for 5 minutes. Reassemble all parts of the nebulizer and connect

it to the air outlet (Pic. G) and switch on the device for 10-15 minutes.

 Do not boil or autoclave the air hose the masks.

## **NEBULIZER REPLACEMENT**

The nebulizer must be replaced after a long period of inactivity, in case it shows deformations or breakings, or when the nebulizer nozzle is obstructed by dry medicine, dust, etc.

DeVilbiss recommends to replace the nebulizer after a period between 6 months and 1 year, depending on the usage.

**Use original nebulizers only.**

## **FILTER REPLACEMENT**

In normal usage conditions, the air filter must be replaced approximately after 100 working hours or once a year. DeVilbiss recommends to periodically check the air filter (10 – 12 treatments) and, if the filter shows a grey or brown colour or is wet, replace it. Remove the filter (Pic. A-4) and re- place it with a new one. Do not attempt to clean the filter to reuse it.

 **Only use original filters.**  
**Do not use the device without the filter.**

## TROUBLESHOOTING

Problem	Possible solution
The device fails to switch on	<ul style="list-style-type: none"> <li>• Make sure the battery is charged (if empty, charge it as indicated in this owner's manual and switch the device on after 5 minutes.</li> <li>• Make sure the plug is properly inserted in the power outlet.</li> </ul>
The device does not nebulize or nebulizes weakly.	<p>Make sure the air hose fittings are properly inserted in the compressor and the nebulizer connections.</p> <p>Make sure the nebulizer is not empty or that it has been filled with the correct amount of medication (MAX. 6 ml). Make sure that the nebulizer nozzle is not clogged.</p>
The device is noisy.	Make sure that the filter is mounted properly.

In the event that the equipment does not resume correctly, please contact DeVilbiss customer service.

## MAINTENANCE AND REPAIRS

In case of failure, contact **qualified personnel authorised** by DeVilbiss. **Never open the device.** The device has no user-serviceable parts within and does not need internal maintenance or lubrication.

## TECHNICAL SPECIFICATIONS

**Model:** MQ6010 AirForce Mini

**Electrical power supply:** Internal lithium polymer (Li-Pol) battery 3.7 V-2570 mAh

**Power supply adapter:**

**Input:** 100-240 V AC, 50-60 Hz 0.5 A

**Output:** 5V DC 2A

**Minimum filling volume:** 2 ml

**Maximum filling volume:** 6 ml

**Nebulization rate:** 0.25 ml/min

**Operating Pressure:** 0.25 ~ 0.5 bar

**Weight:** 217 g

**Dimensions:** 45 x 68 x 111 mm

**Noise level:** 45 dBA

**Expected service life of the device:** 400 hours.

- Class II device referred to protection against electric shock.
- Nebulizer, mouthpiece, masks, and fork are type BF applied parts
- Device for intermittent use (30 min on - 30 min off)
- Device not suitable for anaesthesia and lung ventilation systems.
- Device not suitable for use in the presence of any anaesthetic mixture inflammable with oxygen or nitrogen peroxide.

*Technical specifications subject to change without notice.*

*The internal batteries cannot be replaced by the user.*


## AEROSOL CHARACTERISTICS IN ACCORDANCE WITH REGULATION EN 13544-1 APPENDIX CC

Aerosol output: 0.15 ml

Particle size (MMAD): 4.12  $\mu\text{m}$

### Standards applied:

Electric safety standards CEI EN 60601-1 Electromagnetic compatibility according to CEI EN 60601-1-2 The device is a Class IIa medical device according to Directive 93/42/EEC “Medical devices”.

## ENVIRONMENTAL CONDITIONS

### Storage and Transport

Temperature: MIN  $-20^{\circ}\text{C}$  - MAX  $+60^{\circ}\text{C}$

Humidity: MIN 10%RH - MAX 95%RH

Atmospheric pressure: 700 hPa - 1060 hPa

### Operation

Temperature: MIN  $+10^{\circ}\text{C}$  - MAX  $+40^{\circ}\text{C}$

Humidity: MIN 10%RH - MAX 95%RH

Atmospheric pressure: 700 hPa - 1060 hPa

## SYMBOLS AFFIXED TO THE DEVICE

	ON
	OFF
	TYPE BF.
	CLASS II
	ALTERNATING CURRENT
	ATTENTION
	PRODUCT CODE
	SERIAL NUMBER
	NUMBER THAT IDENTIFIES THE PRODUCTION BATCH
	MANUFACTURER
	TEMPERATURE
	READ THE INSTRUCTIONS
<b>CE 0123</b>	European Community Directive 93/42/EEC
<b>IP22</b>	Protection rating of device casing against the penetration of external solid and liquid agents

**CORRECT DISPOSAL OF THE PRODUCT** (waste electrical & electronic equipment)  
(Applicable in the European Union and countries with separate collection systems).

The symbol on the product/ documents indicates that the product complies with legislation on electrical and electronic equipment and it must not be disposed of as household waste. The user is responsible for taking the device to an appropriate waste disposal centre at the end of its service life, subject to administrative penalties and fines, provided for by the current legislation on waste. For more detailed information on available collection systems, refer to the local waste disposal service.

The product contains a rechargeable battery covered by the European Directive 2006/66/EC,

which cannot be disposed of as regular municipal waste. Please inform yourself about the local rules on separate collection of batteries.

## **TERMS AND CONDITIONS OF WARRANTY**

The device is guaranteed for 2 years from the date of original purchase against any defect in materials or workmanship.

The warranty consists in the replacement and/or repair, free of charge, of originally defective components.

The warranty does not cover the accessories supplied and the parts subject to normal wear and tear. The device may only be repaired by authorised technical service centres.

The appliance must be sent to DeVilbiss Customer Service for repairs.

The transport costs shall be borne by the user.

Any repair out of warranty shall be borne by the user. The warranty lapses if the device has been tampered with, if the defect was caused by improper use or in case the damage is not due to the manufacturer (accidental fall, incorrect transport etc.).

The warranty does not involve any direct or indirect damages of any kind to people or property during the period of inefficiency of the product.

The warranty is valid from the date of purchase certified by the receipt or invoice that must be attached to the warranty coupon.

If the user fails to submit the warranty coupon duly compiled and validated by the purchase document, he will have no right to avail of service under warranty.

**COUPON TO BE RETURNED IN CASE OF REPAIRS**  
**THE WARRANTY IS ONLY VALID IF YOU ATTACH YOUR RECEIPT**

**Device Type:** \_\_\_\_\_

**Model:** \_\_\_\_\_

**Serial Number:** \_\_\_\_\_

**Date of purchase:** \_\_\_\_\_

**BUYER'S DATA**

**Full Name:** \_\_\_\_\_

**Address:** \_\_\_\_\_

**Telephone:** \_\_\_\_\_

**Description of fault:** \_\_\_\_\_

**Signature for acceptance of the warranty conditions**

\_\_\_\_\_

I authorize the use of the above information according to law 675/96 on Privacy.

### Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment

<b>Emission test</b> Emissions RF CISPR 11	<b>Compliance</b> Group 1	<b>Electromagnetic environment – guidance</b> The AirForce Mini uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions  IEC 61000-3-3	Complies	

### Guidance and manufacturer's declaration - electromagnetic emissions

The AirForce Mini is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment

<b>Immunity test</b>  Electrostatic discharges  IEC 61000-4-2	<b>IEC 60601 Test level</b>  ±8 kV contact  ±15 kV air	<b>Compliance level</b>  ±8 kV contact  ±15 kV air	<b>Electromagnetic environment - guidance</b> Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
--	--	--	---

NOTE: UT is the a.c. mains voltage prior to application of the test level.

### Guidance and manufacturer's declaration - electromagnetic emissions

The AirForce Mini is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrical fast transient/burst IEC 61000-4-4	Power supply lines : $\pm 2$ kV input/output lines: $\pm 1$ kV	Power supply lines : $\pm 2$ kV input/output lines: $\pm 1$ kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): $\pm 1$ kV line(s) to earth: $\pm 2$ kV 100 kHz repetition frequency	line(s) to line(s): $\pm 1$ kV to line(s) to earth: $\pm 2$ kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	0% 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°  0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m 50Hz/60Hz	3 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic emissions

The AirForce Mini is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that it is used in such an environment

Immunity test	IEC 60601	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC61000-4-6	<b>Test level</b> 150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of this device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  <b>Recommended separation distances</b> $d = 0.35 \sqrt{P}$ $d = 1.2 \sqrt{P}$
Radiated RF IEC61000-4-3	10V/m 80% Am at 1kHz	10V/m 80% Am at 1kHz	80 MHz to 800MHz/ $d = 1.2 \sqrt{P}$ 800 MHz to 2.7 GHz/ $d = 2.3 \sqrt{P}$ where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. <sup>a</sup> Interference may occur in the vicinity of equipment marked with the following symbol:  

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in

the location in which this device is used exceeds the applicable RF compliance level above, the Air- Force Mini should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device

Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the AirForce Mini

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the AirForce Mini as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz (out ISM and amateur radio bands) $d=1.2\sqrt{p}$	150 kHz to 80 MHz (in ISM and amateur radio bands) $d=0.6\sqrt{p}$	80MHz to 800MHz $d=1.2\sqrt{p}$	800MHz to 2.7GHz $d=2.3\sqrt{p}$
0,01	0.12	0.06	0.12	0.23
0,1	0.38	0.19	0.38	0.73
1	1.2	0.6	1.2	2.3
10	3.8	1.9	3.8	7.3
100	12	6	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## Guidance and manufacturer's declaration - electromagnetic emissions

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the AirForce One should assure that it is used in such an environment

Radiated RF  
IEC61000-4-3  
(Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipments

Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 –390	TETRA 400	Pulse modulation b) 18 Hz	1,8	0,3	27
450	380 –390	GMRS 460, FRS 460	FM c) $\pm 5$ kHz deviation 1 kHz sine	2	0,3	28
710	704 – 787	LTE Band 13,17	Pulse modulation b) 217 Hz	0,2	0,3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18 Hz	2	0,3	28
870						
930						
1720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0,3	28
1845						
1970						
2450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0,3	28
5240	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0,2	0,3	9
5240						
5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
- c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E=6/d \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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Globalcare Medical Technology Co., Ltd
Donawa Lifescience Consulting Srl
DeVilbiss Healthcare GmbH