**drive**

**DeVilbiss®**  
HEALTHCARE

## Drive DeVilbiss® Vacu-Aide® Suction Unit

CE 0044

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## SYMBOL DEFINITIONS

	Attention, consult instruction guide		Center positive polarity indicator		Battery charging
	It is mandatory to read Operating Instructions prior to use. This symbol has a blue background on the product label.		Type BF equipment-applied part		Low battery
	Date of manufacture		"On" compressor		Keep dry
	Manufacturer		"Off" compressor (external battery charging)		Do not get wet
	Direct current		External power		European Representative
	Alternating current		Catalog Number		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
	Container tubing vacuum connector		Choking Hazard – Small parts not for children under 3 years or any individuals who have a tendency to place inedible object in their mouths.		
<b>IP22</b>	Protected against solid material of $\geq 12.5$ mm and vertically falling water drops when enclosure is tilted up to 15°				
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment (WEEE)				

## IMPORTANT SAFEGUARDS

When using electrical products, especially when children or pets are present, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

**DANGER**– Urgent safety information for hazards that will cause serious injury or death.

**WARNING**– Important safety information for hazards that might cause serious injury.

**CAUTION**– Information for preventing damage to the product.

**NOTE**– Information to which you should pay special attention.

## READ ALL INSTRUCTIONS BEFORE USING THIS DEVICE.

THE USER SHALL BE PROPERLY TRAINED PRIOR TO OPERATING THE DEVICE.

## SAVE THESE INSTRUCTIONS.

### DANGER

To reduce the risk of electrocution:

1. Do not use while bathing.
2. Do not place or store product where it can fall or be pulled into a tub or sink.
3. Do not place in or drop into water or other liquid.
4. Do not reach for a product that has fallen into water. Unplug immediately.

### WARNING

To reduce the risk of burns, electrocution, allergy, fire or injury to persons:

1. Close supervision is necessary when this product is used by, on, or near children or physically incapacitated individuals.
2. Use this product only for its intended use as described in this guide.
3. Keep the power cord away from heated surfaces.
4. Never use while drowsy or asleep.

5. Do not cover the unit or the AC to DC adapter while power is applied.
6. Never operate this product if
  - a. It has a damaged power cord or plug.
  - b. It is not working properly.
  - c. It has been dropped or damaged.
  - d. It has been dropped into water.
 Instead return the product to an authorized DeVilbiss Healthcare service center for examination and repair.
7. Check the collection container for cracks before each use. Do not use if cracked.
8. Universal precautions must be observed while providing care or handling equipment for persons with an infectious condition.
9. Do not place the suction unit in a position that would make it difficult to disconnect the AC power cord.
10. When the device is used under extreme operating conditions and in the carrying case, the surface temperature of the device may exceed 73°C.
11. The suction tube shall be kept isolated from other parts of the system to confirm that it does not exceed 43°C prior to patient contact.
12. The catheter to be used shall be carefully inserted to prevent strangulation/choking.
13. Don't let children play with the tubing as this may lead to strangulation.
14. Consult the physician in case of any allergic reaction due the use of the device.
15. Keep the suction pump out of the reach of children and pets and away from pests.

**CAUTION**– Do not connect 6' (1.8m) patient tubing directly to suction unit; damage to suction unit will occur.

**NOTE**– All the materials used in the suction pump are Latex free.

## INTERNATIONAL TRAVEL

The 7314 series is equipped with an AC to DC adapter allowing operation on any AC voltage (100-240 VAC, 50/60 Hz). However the correct power cord must be used to connect to adaptable wall power. AC power is removed by disconnecting the power cord from the AC outlet receptacle.

**NOTE**– Check power cord for adaptability before using.

## INTRODUCTION

Your DeVilbiss suction unit is a compact medical suctioning device which has been designed for reliable, portable operation. Following the recommended operating and maintenance procedures outlined in this instruction guide will maximize the life of this product. This guide provides an overview of the suction unit and operation. Your equipment provider should give more detailed instructions related to patient care.

### Intended Use Statement

The device is to be used to remove fluids from the airway or respiratory support system and infectious materials from wounds. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection container. The fluids are trapped in the collection container for proper disposal. It is for use on the order of a physician only.

### Contraindications (specific situation in which the device shall not be used)

The Vacu-Aide QSU should not be used for:

- thoracic drainage
- nasogastric suction

### DANGER

**The DeVilbiss Vacu-Aide is a vacuum suction device designed for the collection of nonflammable fluid materials in medical applications only. Improper use during medical applications can cause injury or death. For all medical applications:**

1. All suctioning should be done in strict accordance with appropriate procedures that have been established by a licensed medical authority.
2. Some attachments or accessories may not fit the tubing supplied. All attachments or accessories should be checked prior to use to assure proper fit.

## ACCESSORY/REPLACEMENT ITEMS

The following items can be purchased separately as accessories or replacement items for your 7314 Series DeVilbiss Suction Unit:

Description	Part No.	Description	Part No.
6' Patient Tubing (USA)	SUCP TUBING 72	AC to DC adapter/charger	7314P-613
6' Patient Tubing (International)	6305D-611	Power cord for USA	DV51D-606
Collection Container Kit (Internal filter cartridge, splash guard, 800 ml container, 4-1/2" (11.43 cm) and 6' tubing package)	7305D-633	Power cord for Continental Europe	DV51D-607
800 ml disposable container with internal filter cartridge, splash guard and 4-1/2" (11.43 cm) tubing (48 each)	7305D-632	Power cord for UK	DV51D-608
Filter cartridge (12 pack) (For Disposable Container)	7305D-635	Power cord for Australia	DV51D-609
Collection Container Kit (1200 ml reusable container, external bacteria filter, elbow, 4-1/2" (11.43 cm) tubing)	7314D-603	Power cord for Brazil	DV51D-612
1200 ml reusable container (external bacteria filter, elbow, 4-1/2" (11.43 cm) tubing) (6 pack)	7314D-604	Power cord for Japan	DV51D-613
External bacteria filter (non-sterile) (12 pack) For reusable container	7305D-608	Power cord for China	DV51D-614
Carrying case	7314D-606	Power cord for Argentina	180-0006-011
12V DC power cord (1 each)	7304D-619		


**NOTE**– The use of electrical cables and accessories other than those specified in this manual or referenced documents may result in increased electromagnetic emissions from the product or decreased electromagnetic immunity of the product.

## IMPORTANT PARTS

### 7314 Series DeVilbiss Vacu-Aide QSU Suction Unit


#### with Disposable Container and Filter Cartridge

1. Vacuum gauge
  2. Vacuum regulator knob
  3. 6' patient tubing
  4. Patient tubing connector
  5. Disposable container with lid (float shut off incorporated into lid) and filter cartridge
  6. Filter cartridge with 4-1/2" (11.43 cm) tubing
  7. DC power input (on side)
  8. Power switch
  9. LED power lights
- AC to DC adapter (not shown)  
DC power cord (not shown) optional  
Internal rechargeable battery (not shown) 7314P series only  
Carrying case (not shown) 7314P series only

#### with Reusable Container and External Bacteria Filter

1. Vacuum gauge
  2. Vacuum regulator knob
  3. 4-1/2" (11.43 cm) connection tubing
  4. External bacteria filter
  5. Patient tubing connector
  6. 6' patient tubing
  7. Lid
  8. DC power input (on side)
  9. Power switch
  10. LED power lights
- AC to DC adapter (not shown)  
DC power cord (not shown) optional  
Internal rechargeable battery (not shown) 7314P series only  
Carrying case (not shown) 7314P series only

#### Disposable Container with Filter Cartridge and Splash Guard

1. 4-1/2" (11.43 cm) connection tubing
2. Filter cartridge (Do not get wet)
3. Lid
4. Jar
5. Splash guard
6. Patient tubing connector


#### Reusable Container with External Bacteria Filter

1. 4-1/2" (11.43 cm) connection tubing
2. Lid with o-ring
3. Overflow valve
4. Jar
5. Patient tubing connector
6. Connection elbow
7. Bacteria filter


## SET-UP & OPERATION

1 Fully charge battery for 17 HOURS. (7314P series only)

2 If applicable, ensure splash guard is securely attached to inside of lid over filter cartridge.

3 Securely attach lid to container.

4 Insert container into holder and gently push into place.  
**NOTE** – Do not use excessive force. Pushing container down too hard could cause potential leak or loss of suction.

5A Disposable Container w/ internal filter cartridge  
**Connection:** Attach 4-1/2" (11.43 cm) tubing from filter cartridge to tubing connector on unit.

5B Reusable Container w/external Bacteria Filter Connection:  
Connect either end of the 4-1/2" (11.43 cm) tubing to the tubing connector then connect the other end to the bacteria filter. Ensure that the clear side of the bacteria filter is toward elbow and bottle when installing. Do not reverse direction of filter. The bacteria filter should then be connected to the 90° elbow connection, and the elbow should be connected to the top of the container lid where it says <Vacuum>.

6 Attach 6' patient tubing to container lid at outlet labeled <Patient>.

7 Ensure power switch is "off".

**CAUTION**– Do not connect 6' (1.8m) patient tubing directly to suction unit; damage to suction unit will occur.

**NOTE**–Inspect suction tubing and container for leaks, cracks, etc. and assure that all connections are secure and without leaks before using.

7A 7314P - Select desired power source. (Skip step 8 if using internal battery power.)


7B 7314D series (non-battery label)


**NOTE**– The 7314D series is not factory equipped with an internal rechargeable battery. 7314P series is factory equipped with an internal rechargeable battery and all information regarding battery operation in this guide is applicable.

8 If using AC or DC power, plug the small connector into the DC power input on the side of unit.

9 Plug the other end into an AC wall outlet or DC receptacle.  
**NOTE**– The AC adapter may become warm to the touch during charging or running of the unit. This is normal.

10 Turn the unit "on".

11 Adjust the suction level.

12 Verify suction level.  
**NOTE** - Always verify suction level before beginning by occluding open end of patient tubing while observing gauge. Adjust knob to desired level.

### WARNING

If the unit is not receiving power from an external source or the battery was not recharged, the low battery indicator light will remain on and the performance of the unit will drop off rapidly. Switch to another power source immediately after the low battery light appears to avoid an interrupted suction procedure.

**NOTE**– The vacuum (negative) pressure should be set as directed by a physician or other appropriate health care professional. Published expert consensus suggests the suction vacuum pressure should be set as low as possible to effectively clear secretions.<sup>1</sup> For tracheal and/or endotracheal suctioning, a vacuum (negative) pressure of less than 100 mmHg in infants, less than 120 mmHg in children and less than 150 mmHg in adults has been recommended.<sup>2</sup>

**NOTE**– Gauge is for reference only. If the unit sustains a severe drop, accuracy of the gauge must be checked.

**CAUTION**– When automatic float shut-off is activated, contents of the collection container should be emptied. Further suctioning could cause damage to the vacuum pump.

**CAUTION**– Should fluid be aspirated back into the unit, equipment provider servicing is necessary as possible vacuum pump damage may result.

### WARNING

This device is rated for intermittent operation use as noted in the specifications. Do not operate device with no flow (or suction port blocked). Continuous operation under this condition may cause the unit to thermal shut off or damage the device.


<sup>1</sup> AARC Clinical Practice Guideline. Endotracheal Suctioning of Mechanically Ventilated Patients With Artificial Airways-2020. Respir Care 2010;55(6): 758- 764.

<sup>2</sup> AARC Clinical Practice Guideline. Nasotracheal Suctioning-2004 Revision & Update. Respir Care 2004;49(9): 1080-1084

# BATTERY CHARGING & FILTER MAINTENANCE

## Battery Charging (7314P Series Only)


On 7314P series, the units are equipped with a factory-installed rechargeable battery. The unit will have a light for low battery and charge indication.


1 Ensure power switch is "off".


2 Plug the small connector of the AC or DC adapter into the DC power input.

3 Plug the other end into an AC wall outlet or DC receptacle.

4 Battery charging begins; 17 hours for full charge.

5 Battery charging complete.

### LED Explanations:

**Green**— Illuminated when external power is supplied to unit from an AC or DC power source.

**Yellow**— Battery is being charged. Light will go out when battery is fully charged.

**Red**— Low battery. Seek another power source and charge battery as soon as possible when light illuminates.

**NOTE** – Charge battery for a minimum of 17 hours before first use.

**NOTE** – Fully recharge battery after each use. The unit will continue to float charge the battery after the charge indicator turns off, so leave the unit connected to AC when not in use.

**CAUTION**– Discharging the battery completely will shorten the life of the battery. Do not operate the unit more than a few minutes if the low battery indicator light is lit. Recharge as soon as possible.

**NOTE** – Unit run time will decrease as the battery ages.

**NOTE** – Unit run time will also be reduced by letting the battery sit at a discharged state for extended periods.

**STORAGE NOTE** - Battery should be charged for a minimum of 17 hours prior to storage, and at least once every 6 months. Important - If battery recharge is delayed beyond 6 months, battery may be able to provide full run time after completion of 3 full charge and discharge cycles.

**NOTE**– A fully charged battery will provide approximately 60 minutes of continuous operation at a zero vacuum level (free flow). Operation time will decrease with higher vacuum levels.

**NOTE**– When charging the battery, use an external power source and verify that the charge light illuminates when the unit is in the "Off" position. If the battery does not charge, please be sure the model you are using has a battery installed prior to contacting your authorized DeVilbiss Healthcare provider.


**NOTE**– The internal rechargeable battery is sealed lead-acid. Contact your local authorities for instruction on proper disposal.

**NOTE**– Do not connect the AC adapter to an outlet controlled by a switch to ensure power is supplied to unit at all times.

**NOTE**– Do not connect the DC power cord to an outlet that is not constantly energized.

## Changing Filter Cartridge (single-patient use) Disposable Container

Change filter cartridge if overflow occurs or every two months, whichever comes first.

1 Turn unit "off".

2 Remove filter cartridge and 4-1/2" (11.43 cm) tubing.

3 Install new cartridge and attach tubing.


**NOTE**– Do not substitute any other material for this filter cartridge. Substitution may lead to contamination or poor performance; use only DeVilbiss filter cartridges.

**NOTE**– The filter cartridge contains a hydrophobic filter. If the filter media becomes wet, air flow will be stopped. The filter cartridge must then be replaced. Do not remove filter media from filter cartridge.


**NOTE**– Filter cartridges are included with each disposable container. They are also available separately (7305D-635 12/pack).

## Changing Bacteria Filter (single-patient use) Reusable Container

Change bacteria filter if overflow occurs or every two months, whichever comes first.

1 Turn unit "off".

2 Remove filter by disconnecting it from suction unit and lid assembly.

Replace with a clean DeVilbiss bacteria filter (non-sterile) and remount to suction unit and lid. Ensure that the clear side of the bacteria filter is toward elbow and bottle when installing. Do not reverse direction of filter. Additional filters (7305D-608 12/pack) may be purchased from your authorized DeVilbiss Healthcare provider.

**NOTE**– Do not substitute any other material for this bacteria filter. Substitution may lead to contamination or poor performance; use only DeVilbiss filters.

**NOTE**– Bacteria filter must be changed between patients.

## Service Life

- Service life of the suction unit is 5 years.
- Service life of the battery is 200 discharge cycles.
- Service life of the disposable container is 30 dishwasher cycles.
- Service life of the reusable container is 30 autoclave cycles at 121°C.

# CLEANING INSTRUCTIONS

## WARNING

To prevent possible risk of infection from contaminated cleaning/disinfection solutions, always prepare fresh solution for each cleaning cycle and discard solution after each use.


**NOTE**– Disinfection information is based on AARC Clinical Practice Guideline Suctioning of the Patient in the Home.

## Disposable Collection Container (with Internal Filter Cartridge)


**NOTE**– The disposable collection container and lid are meant for single-patient use.

Turn unit "off" and allow vacuum to drop.

Disconnect from power source.

Disconnect tubing and remove container from holder.

Carefully remove lid and empty contents.

Remove filter cartridge and 4-1/2" (11.43 cm) tubing and set aside.

Filter **MUST NOT** get wet. The filter material cannot be removed from the elbow.


Wash container, lid and splash guard in warm water/dishwashing solution. Rinse with clean, warm water.


Soak in 1 part vinegar ( $\geq 5\%$  acetic acid concentration) to 3 parts water (131°F-149°F/55°C-65°C) solution for 60 minutes. Rinse with clean, warm water and air dry.

**WARNING**- Do not remove float ball from lid. If removed, float ball may pose a choking hazard.


**NOTE**– The disassembled container may also be washed in a dishwasher, top shelf only, using a cycle with a water temperature between 131°F-149°F/55°C-65°C.


## Reusable Collection Container (with External Bacteria Filter)


Turn unit "off" and allow vacuum to drop.

Disconnect from power source.

Disconnect tubing and remove container from holder.

Carefully remove lid and empty contents.

Remove bacteria filter, elbow, and 4-1/2" (11.43 cm) tubing and set aside.

Wash jar, lid, o-ring, and overflow valve in a solution of warm water with a mild, liquid detergent (e.g. Dawn or Palmolive) and rinse with clean, warm tap water. Then disinfect using one of the following methods.

### For single patient use:

1. Soak in 1 part vinegar ( $\geq 5\%$  acetic acid concentration) to 3 parts water (131°F-149°F/55°C-65°C) solution for 60 minutes. Rinse with clean, warm water and air dry in a clean environment.
2. Soak with a commercial (bacterial-germicidal) disinfectant. Follow disinfectant manufacturer's recommended dilution rates and instructions carefully.

### For multi patient use:

1. After parts are completely dry, place jar and lid in autoclave with the open end down. Ensure parts are not touching. Run one sterilization steam cycle at 250°F (121°C) for 15 minutes. **NOTE**-Jar is guaranteed up to 30 cycles of autoclave sterilization at the indicated conditions.
2. Dispose of and replace filter, tubing and elbow between patients.

## 6' Patient Tubing (single-patient use)


1 Disconnect from lid.

2 Rinse thoroughly by running warm tap water through it.

3 Follow by soaking in a solution of 1 part vinegar ( $\geq 5\%$  acetic acid concentration) to 3 parts water (131°F-149°F/55°C-65°C) for 60 minutes. Rinse with clean, warm water and air dry.



4 Keep outer surface clean by wiping with clean, damp cloth.

## AC to DC Adapter

1 Disconnect AC to DC adapter from unit and from power source.

2 Wipe AC to DC adapter housing and cords with a dry cloth.


## Suction Unit (single-patient use)

**CAUTION**– Do not submerge in water as this will result in damage to the vacuum pump.

**NOTE**– Do not use any cleaners or disinfectants that contain ammonia, benzene and/or acetone to clean the unit.

1 Turn unit "off" and allow vacuum to drop.

2 Disconnect from power source.

3 Wipe the housing with a clean cloth and any commercial (bacterial-germicidal) disinfectant.

## Carrying Case (single-patient use)

1 Wipe with clean cloth dampened with detergent or disinfectant.

## Suction Unit (multi-patient use)

### Device Cleaning and Disinfection When There is a Patient Change

DeVilbiss Healthcare recommends that at least the following procedures be carried out by the manufacturer or a qualified third party between uses by different patients.

**NOTE**– When the device is used as per instructions under normal conditions the interior of the device is protected from exposure to pathogens by the in-line filter on the collection container, hence no disinfection of internal components is necessary.

**NOTE**– If the device is used without an in-line filter then the interior of the device may be exposed to pathogens and the device cannot be disinfected.

**NOTE**– If the following described complete processing of the device by an appropriately trained individual is not possible, the device should not be used by another patient.

1. Use disinfectants safely. Always read the label and product information before use.
2. Always wear personal protective equipment when performing this procedure. Use suitable gloves and safety glasses. Cover exposed skin on arms to prevent accidental contact with bleach solution that has been applied to the device. Perform disinfection in a well-ventilated area.
3. Dispose of all accessories that are not suitable for reuse. This includes but may not be limited to the collection container, filter, tubing and carrying case.
4. With the power switch in the "Off" position, disconnect the device from all external power sources.
5. Visually inspect the device for any damage, missing parts, etc.
6. Clean the exterior of the device with a clean lint-free cloth. Heavy soil should be removed with a clean lint-free cloth dampened with water. A soft bristled brush dampened with water can be used to remove stubborn soil. Dry the device using a clean lint-free cloth if water was used to remove soil.
7. Use 5.25% chlorine bleach (Clorox Regular Liquid Bleach or equivalent). Mix one (1) part bleach with four (4) parts water in an appropriate clean container. This ratio produces a one (1) part bleach to five (5) total parts solution (1:5). The total volume (amount) of solution required is determined by the number of devices in need of disinfection. **NOTE**– An alternate suitable commercial (bacterial-germicidal) disinfectant that meets the requirements listed in the NOTE below and is used as per the disinfectant manufacturer's recommended dilution rates and instructions may also be used.

**CAUTION**– Do not submerge the device in water as this will result in damage to the vacuum pump.

**NOTE**– Do not use any cleaners or disinfectants that contain ammonia, benzene and/or acetone to clean the device.

8. Apply the bleach solution in an even manner to the enclosure and power cord using a clean lint-free cloth. The cloth should be dampened only and not dripping of solution. Do not use a spray bottle to apply the solution. Do not saturate the device with the solution. Avoid over-saturating the enclosure seams so that no solution residue builds up in these areas.
9. Exposure time of the disinfectant solution should be 10 minutes minimum to 15 minutes maximum.
10. After the recommended exposure time, all exterior surfaces of the device are to be wiped with a clean lint-free cloth dampened with drinking quality water no warmer than room temperature. Dry the device with a dry, clean lint-free cloth. This is to remove residue that may stain or leave a film on the device, especially after repeated disinfections.

## Disinfection

**NOTE** – The disinfection process can only be completed by the manufacturer or by an appropriately trained individual.

	<b>Recommended disinfection interval</b>	<b>Number of disinfection cycles</b>	<b>Compatible disinfection method</b>
Enclosure, power cord	Between patients	20	1:5 chlorine bleach (5.25%) and water solution, Mikrobac forte, Terralin Protect
Collection container, filter, tubing and carrying case	Do not clean, replace between patients	N/A	N/A

## PROVIDER'S NOTES

No routine calibration or service is required provided the device is used in accordance with the manufacturer's directions. In case of a change of patient, the device must be reconditioned to protect the user. Reconditioning must only be carried out by the manufacturer or service provider. Between patients:

1. Visually inspect unit for any damage, missing parts etc.
2. Ensure that unit and accessories are clean.
3. Use an independent vacuum gauge to verify the unit provides the proper vacuum level as stated in Specifications.
4. Discard and replace collection container, filter, tubing, and carrying case between patients.
5. Wipe the surface using a clean cloth dampened with disinfectant.

## TROUBLESHOOTING

### DANGER

**Electric shock hazard. Do not attempt to open or remove cabinet, there are no user-serviceable internal components. If service is required, return unit to a qualified DeVilbiss Healthcare provider or an authorized service center. Opening or tampering with the unit will void warranty.**

**NOTE**– Your DeVilbiss Suction Unit contains no user-serviceable parts. If you believe your unit is not working properly, **BEFORE YOU RETURN IT TO your DeVilbiss Healthcare PROVIDER WHERE YOU PURCHASED IT**, please take a few moments to check for these possible causes:

<b>PROBLEM</b>	<b>ACTION</b>
<b>Unit does not turn on when external power is connected. Green external power light does not illuminate.</b>	<ol style="list-style-type: none"><li>1. Check power sources and connections.</li><li>2. Ensure wall outlet is live by plugging in a lamp.</li></ol>
<b>Pump runs, but there is no suction.</b>	<ol style="list-style-type: none"><li>1. Check that all tubing is connected properly.</li><li>2. Check tubing connections for breaks or leaks.</li><li>3. Ensure that float shut-off in collection container is not activated or filter cartridge occluded.</li><li>4. Check for leaks or cracks in collection container assembly.</li></ol>
<b>Low suction.</b>	<ol style="list-style-type: none"><li>1. Use vacuum regulator knob to increase suction level.</li><li>2. Check system for leaks.</li></ol>
<b>Unit does not turn on (no external power is connected). 7314P Series Only</b>	<ol style="list-style-type: none"><li>1. Check that battery is fully charged and/or charge battery.</li></ol>
<b>Battery will not charge (external power and charge indicator lights should be illuminated during charge mode) 7314P Series Only</b>	<ol style="list-style-type: none"><li>1. Verify that both external power and charge indicator lights illuminate.</li><li>2. Check power sources and connections.</li><li>3. Ensure wall outlet is live by plugging in a lamp.</li></ol>

## SPECIFICATIONS/CLASSIFICATIONS

Size (including container)	8.3 H x 8.0 W x 8.5 D (21.1 cm x 20.3 cm x 21.6 cm) (not including AC to DC adapter)
Weight (including container)	7314P Series - 6.6 lb. (3 kg) (not including AC to DC adapter) 7314D Series - 4.3 lb. (2.0 kg)
Typical Operating Sound Level	55 dBA
Electrical Requirements	100-240V~, 50/60Hz, 1.2A max $\cdot \text{---} \bullet \text{---} +$ ; 12V $\text{---} \text{---} \text{---}$ ; 33 W Max
Vacuum Range	50-550 mm Hg +/- 10%*
Air Flow @ Pump Inlet:	27 LPM (free flow) typical (may be less when running from internal battery)* ( <b>NOTE</b> – Using Flovac Liners may impair the device's performance.)
Disposable Collection Container Capacity	800 ml (cc)
Reusable Collection Container Capacity	1200 ml (cc)
<b>Warranty</b>	Two-years limited, excluding internal battery (7314P series only) and collection container
Internal Battery (7314P Series Only)	90-day
<b>Approvals</b>	IEC 60601-1:2005+AMD1:2012; EN 60601-1-2: 2015; ISO 10079-1:2015; IEC 60601-1-2:2014; IEC 60529:1989+AMD1:1999+AMD2:2013; IEC 60601-1-6:2010; IEC 60601-1-11:2015; IP22 CONFORMS TO AAMI STD. ES 60601-1, AAMI STD. HA 60601-1-11, IEC STD. 60601-1-6, ISO STD. 10079-1 CERTIFIED TO CSA STD. C22.2 NO. 60601-1, NO. 60601-1-11, IEC STD. 60601-1-6, ISO STD. 10079-1 7314 Series complies with IEC 60601-1 3rd Edition
<b>Environmental Conditions</b>	
Operating Temperature Range	32°F (0°C) - 104°F (40°C)
Operating Relative Humidity	0-95%
Operating Atmospheric Pressure	10.2 psi (70 kPa) - 15.4 psi (106 kPa)
Storage & Transport Temperature Range	-40°F (-40°C) - 158°F (70°C)
Storage & Transport Relative Humidity	0-95%
Storage & Transport Atmospheric Pressure	7.3 psi (50 kPa) - 15.4 psi (106 kPa)
<b>Equipment Classifications</b>	
With respect to protection from electric shock	Class II and internally powered
Degree of protection against electric shock	Type BF Applied Parts
Degree of protection against ingress of liquids	IP22 and ordinary power supply
Mode of Operation	Intermittent Operation: 30 minutes on, 30 minutes off
Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.	
<b>ISO Classification</b>	
High Flow/High Vacuum	
7314 series - Electrically powered medical suction equipment for non-transport use according to EN ISO 10079-1:2015	

\* Conditions may vary based on altitude above sea level, barometric pressure, and temperature.

## TWO-YEAR LIMITED WARRANTY

The compressor portion of the DeVilbiss Vacu-Aide QSU Suction Unit 7314 Series (excluding internal rechargeable batteries and collection containers) is warranted to be free from defective workmanship and materials for a period of two years from date of purchase. Internal rechargeable batteries are warranted for 90 days. Any defective part(s) will be repaired or replaced at DeVilbiss Healthcare's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your DeVilbiss Healthcare Provider or DeVilbiss Service Department for instructions: U.S.A. 800-338-1988 or 814-443-4881, Europe +49-(0) 621-178-98-0.

**NOTE**– Be sure to retain a dated proof of purchase document to verify unit is within 2-year warranty period.

**NOTE**– This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

### Manufacturer's Note

Thank you for choosing a DeVilbiss Suction Unit. We want you to be a satisfied customer. If you have any questions or comments, please send them to our address on the back cover.

For Service Call Your Authorized DeVilbiss Healthcare Provider:

Phone	Purchase Date	Serial #
-------	---------------	----------

## **ELECTROMAGNETIC COMPATIBILITY INFORMATION**

---

### **WARNING**

 MR Unsafe

Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the Vacu-Aide or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.

Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the Vacu-Aide. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.

This device is suitable for use in home and healthcare environments except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of Electromagnetic DISTURBANCES is high.

### **Cables and maximum lengths of cables**

DC power cord (cigarette lighter adapter) #7304D-619 maximum length = 2.2 meter (7.2 ft)

### **WARNING**

Use of accessories and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

### **WARNING**

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

### **WARNING**

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Vacu-Aide, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result

**DeVilbiss Healthcare LLC**

100 DeVilbiss Drive  
Somerset, PA 15501-2125  
USA  
800-338-1988 • 814-443-4881

**Drive DeVilbiss Healthcare Ltd.**

Sidhil Business Park  
Holmfield, Halifax  
West Yorkshire HX2 9TN  
ENGLAND  
+44 (0) 845 0600 333

**Australian Sponsor:**

**Drive DeVilbiss Australia Pty Limited**

Building F, 2 Hudson Ave  
Castle Hill, NSW 2154  
AUSTRALIA  
+61 02 9899 3144

**Drive DeVilbiss Healthcare France**

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DeVilbiss® and Vacu-Aide® QSU are registered trademarks of DeVilbiss Healthcare.

© 2023 DeVilbiss Healthcare LLC. 01.23 All Rights Reserved.

SE-7314-3 Rev. I

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**DeVilbiss Healthcare GmbH**
DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA