EN **DeVilbiss® Vacu-Aide®
Compact Suction Unit
Instruction Guide**

REF **7310 Series**

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



IEC SYMBOLS

	Attention, consult instruction guide		Center positive polarity indicator		Do not get wet
	It is mandatory to read and understand the operating instructions prior to use. This symbol has a blue background on the product label.		Standby - On/Off		Keep dry
	Direct current		External Power		Catalog Number
	Alternating current		Battery Charging		Date of manufacture
	Type BF equipment-applied part		Low Battery		Manufacturer
IP22	Protected against solid material of ≥ 12.5 mm and vertically falling water drops when enclosure is tilted up to 15° .				European Representative
	Choking Hazard – Small parts not for children under 3 years or any individuals who have a tendency to place inedible objects in their mouth.				
	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.				
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/ EU-Waste Electrical & Electronic Equipment				

IMPORTANT SAFEGUARDS

When using electrical products, especially when children or pets are present, basic safety precautions should always be followed. Read all instructions before using. Important information is highlighted by these terms:

DANGER—Urgent safety information for hazards that will cause serious injury or death.

WARNING—Important safety information for hazards that might cause serious injury.

CAUTION—Information for preventing damage to the product.

NOTE—Information to which you should pay special attention.

READ ALL INSTRUCTIONS BEFORE USING.

THE USER SHALL BE PROPERLY TRAINED PRIOR TO OPERATING THE DEVICE.

SAVE THESE INSTRUCTIONS

DANGER

To reduce the risk of electrocution:

1. Do not use while bathing.
2. Do not place or store product where it can fall or be pulled into a tub or sink.
3. Do not reach in or drop into water or other liquid.
4. Do not reach for a product that has fallen into water. Unplug immediately.

WARNING

To reduce the risk of burns, electrocution, allergy, fire or injury to persons:

1. Close supervision is necessary when this product is used by, on, or near children or physically incapacitated individuals.
2. Use this product only for its intended use as described in this guide.
3. Never operate this product if:
 - a. It has a damaged power cord or plug.
 - b. It is not working properly.
 - c. It has been dropped or damaged.
 - d. It has been dropped into water. Return the product to an authorized DeVilbiss service center for examination and repair.
4. Keep the power cord away from heated surfaces.
5. Never use while drowsy or asleep.
6. Do not cover the unit or the AC to DC adapter while power is applied.
7. Check the collection container for cracks before each use. Do not use if cracked.
8. Universal precautions must be observed while providing care or handling equipment for persons with an infectious condition.
9. Do not place the suction unit in a position that would make it difficult to disconnect the AC power cord.
10. When the device is used under extreme operating conditions and in the carrying case, the surface temperature of the device may exceed 73°C.
11. The suction tube shall be kept isolated from other parts of the system to confirm that it does not exceed 43°C prior to patient contact.
12. The catheter to be used shall be carefully inserted to prevent strangulation/choking.
13. Don't let children play with the tubing as this may lead to strangulation.
14. Consult the physician in case of any allergic reaction due the use of the device.
15. Keep the suction pump out of the reach of children and pets and away from pests.

NOTE– All the materials used in the suction pump are Latex free.

INTERNATIONAL TRAVEL

The suction unit is equipped with an AC to DC Adapter/Charger allowing operation on any AC voltage (100-240 VAC, 50/60 Hz). However the correct power cord must be used to connect to adaptable wall power. AC power is removed by disconnecting the power cord from the AC outlet receptacle.

NOTE– Check power cord for adaptability before using.

INTRODUCTION

Your DeVilbiss Suction Unit is a compact medical suctioning device which has been designed for reliable, portable operation. Because of the small size, light weight and DC operation, the DeVilbiss Suction Unit is ideal for providing suction in the home, in transport with optional DC cord, or if your model has an internal rechargeable battery, the unit can be operated anywhere. Following the recommended operating and maintenance procedures outlined in this instruction guide will maximize the life of this product. This guide provides an overview of the suction unit and operation. Your equipment provider should give more detailed instructions related to patient care.

Intended Use Statement

The device is to be used to remove fluids from the airway or respiratory support system and infectious materials from wounds. The device creates a negative pressure (vacuum) that draws fluids through disposable tubing that is connected to a collection container. The fluids are trapped in the collection container for proper disposal. It is for use on the order of a physician only.

Contraindications (specific situation in which the device shall not be used)

The Vacu-Aide Compact Suction Unit should not be used for:

- thoracic drainage
- nasogastric suction

DANGER


The DeVilbiss Vacu-Aide is a vacuum suction device designed for the collection of nonflammable fluid materials in medical applications only. Improper use during medical applications can cause injury or death. For all medical applications:

1. All suctioning should be done in strict accordance with appropriate procedures that have been established by a licensed medical authority.
2. Some attachments or accessories may not fit the tubing supplied. All attachments or accessories should be checked prior to use to assure proper fit.

IMPORTANT PARTS

Inspect the suction unit and all parts before use.

1. Vacuum Regulator Knob (on side)
2. Bacteria Filter (non-sterile)
3. Vacuum Inlet Port (behind filter)
4. Display Panel (top of unit)
5. 725 ml Reusable Container with Lid/Elbow Assembly
6. DC Power Input (on side)
7. Battery Door
8. Unit Carry Handle/Catheter Holder
9. 6' (1.8m) Patient Tubing (not shown)
10. AC to DC Adapter/Charger (not shown)
11. 12V DC Power Cord (not shown)
12. High Capacity Rechargeable Battery (not shown)
13. Carry Case with Shoulder Strap (not shown)

ACCESSORY/REPLACEMENT ITEMS

The following items can be purchased separately as accessories or replacement for your suction unit.

Description	Part No.	Description	Part No.
External bacteria filter (non-sterile) (12 pack) For reusable container and with external filter	7305D-608	Line Cord for USA	DV51D-606
Battery Door	7310P-602	Line Cord for EU	DV51D-607
AC to DC Adapter/Charger	7314P-613	Line Cord for UK	DV51D-608
6' Patient Tubing (USA)	SUCP TUBING 72	Line Cord for Australia	DV51D-609
6' Patient Tubing (International)	6305D-611	Line Cord for Japan	DV51D-613
Carrying Case w/Shoulder Strap	7310P-606	Line Cord for Brazil	DV51D-612
725 ml Reusable Container Kit (Jar, Lid/Elbow Assembly, Filter)	7310P-603	Line Cord for China	DV51D-614
12V DC Power Cord	7304D-619	Line Cord for Argentina	180-0006-011

NOTE– The use of electrical cables and accessories other than those specified in this manual or referenced documents may result in increased electromagnetic emissions from the product or decreased electromagnetic immunity of the product.

SET-UP

1. Your unit was shipped with the battery disconnected. Connect battery and charge before first use. Refer to Battery Replacement section for instructions.
2. Place clear side of filter labeled IN directly into container's lid.
3. With filter attached to container, directly insert into the front port. Both the filter and the bottom support secure the container to the unit (FIG. 1 & 2).
4. The 6' (1.8m) patient tubing can now be connected to the container lid at the outlet labeled <Patient>.
5. Ensure that all connections are secure before use to prevent leaks.
6. Verify that the unit is at the desired suction level before beginning patient suctioning.

HOW TO OPERATE YOUR SUCTION UNIT

Control Panel Symbols

	On/Off
	External power: Supplied from AC or DC power cord. Illuminates in GREEN when power is supplied.
	Battery charging: Illuminates in YELLOW. The light will go out when the battery is fully charged.
	Low battery: illuminates in RED when battery reaches a discharge state.
mmHg	Suction level setting: Scale/strength illuminates in GREEN. This scale shows the level of suction strength in mmHg.

Power Source Options

AC OPERATION - Plug the 90° power connector of AC to DC adapter/charger into DC power input (FIG 3) and attach the line cord. Plug the AC end into a wall-outlet power source.

12V DC OPERATION - Plug the small 90° power connector of DC power cord into DC power input (FIG 3). Plug large end of cord into 12V DC power receptacle of vehicle.

BATTERY OPERATION - Your unit is equipped with a high capacity rechargeable battery. For initial charge on new unit, fully charge the battery for a minimum of 5 hours (see Battery Charging). To operate unit from the rechargeable battery, ensure that no external power sources are plugged into the suction unit.

NOTE– *During charging or operating, the power supply may become warm to touch; this is normal.*

WARNING

If the unit does not receive external power or the battery has not been charged, the low battery indicator light will remain on and the performance of the unit will drop rapidly. If Low Battery symbol lights, switch to another power source to avoid an interrupted suction procedure.

How to adjust the vacuum level

1. Once power source is selected, turn the unit on by pressing the “On” button. The GREEN light, indicating external power, will remain lit when external power is connected.
2. Occlude (block) the patient end of the tubing, then adjust vacuum level from 50-500 mmHG +/- 10% by turning the vacuum regulator knob clockwise to increase and counter-clockwise to decrease the vacuum (FIG 4). Release and occlude once more to confirm setting. The desired level of vacuum can be viewed on the LED display.

NOTE– *The LEDs have two brightness levels. As the vacuum level is adjusted, the LEDs will illuminate in progression. When an LED is at half brightness, it indicates that the vacuum level is halfway between the previous fully lit LED and the half brightness LED. EXAMPLE: If the 150 mmHg LED is fully illuminated and the 200 mmHg LED is at half brightness, the suction level is 175 mmHg. When the 200 mmHg LED illuminates at full brightness, the unit has reached 200 mmHg.*

3. Connect suction tip or catheter as appropriate.

NOTE– *The vacuum (negative) pressure should be set as directed by a physician or other appropriate health care professional. Published expert consensus suggests the suction vacuum pressure should be set as low as possible to effectively clear secretions.¹ For tracheal and/or endotracheal suctioning, a vacuum (negative) pressure of less than 100 mmHg in infants, less than 120 mmHg in children and less than 150 mmHg in adults has been recommended.²*

NOTE– *If the unit does not maintain vacuum, refer to Troubleshooting.*

¹ AARC Clinical Practice Guideline. Endotracheal Suctioning of Mechanically Ventilated Patients With Artificial Airways-2020. Respir Care 2010;55(6): 758- 764.

² AARC Clinical Practice Guideline. Nasotracheal Suctioning-2004 Revision & Update. Respir Care 2004;49(9): 1080-1084

NOTE– Suction ceases if liquid level reaches float shut-off valve located on underside of container lid.

CAUTION– Further suctioning may cause damage to the vacuum pump and voids warranty. Equipment service is required if fluid content is aspirated back into the unit.

WARNING

This device is rated for intermittent operation use as noted in the specifications. Do not operate device with no flow or suction port blocked. Continuous operation under this condition may cause the unit to thermal shut off or damage the device.

BATTERY CHARGING

The Vacu-Aide Compact series is equipped with a factory-installed high capacity rechargeable battery. Connect battery and charge before first use (refer to Battery Replacement).

1. Connect the unit to either an AC or DC power source.
2. The green external power light should illuminate; the yellow charge indicator will remain lit while the battery is charging.
3. Ensure that the yellow charging light is illuminated when charging begins. As the battery nears a full charge, the yellow battery charging light may flash for several minutes. This is normal. If your unit does not hold a charge, check that the yellow light turns on when external power is applied with the power button "Off." If problems persist, contact your medical equipment provider.

NOTE– A fully discharged battery will require up to 5 hours (depending on depth of discharge) of charging to reach full capacity.

If unit is not in use for extended periods, the battery should be recharged a minimum of every 6 months. A fully charged battery will provide approximately 45-60 minutes of continuous operation at zero vacuum level (free flow).

CAUTION– Completely discharging the battery will shorten the battery life. Do not operate the unit for more than a few minutes if the low battery indicator light is lit. Recharge the battery as soon as possible.

BATTERY REPLACEMENT

1. Using a coin or straight-blade screwdriver rotate latch to unlocked position (FIG. 5).
2. Remove door by pulling up on latch (FIG. 6).
3. Remove battery from compartment and unplug connector from circuit board (FIG. 7 & 8).
4. Install new battery by reversing the above steps.
5. Fully charge battery before using.

CLEANING INSTRUCTIONS

WARNING

To prevent possible risk of infection from contaminated cleaning /disinfection solutions, always prepare fresh solution for each cleaning cycle and discard solution after each use.

Collection Container:

NOTE– Elbow is permanently attached to the lid and cannot be removed.

1. Push power button to turn unit off. Wait for vacuum level to drop.
2. Disconnect external power source from input receptacle on unit (if applicable).
3. Remove container from unit by disconnecting the patient tubing and filter.

NOTE– The reusable collection container and lid are single-patient use and should be replaced every 1 to 2 months after a minimum of 30 cleaning cycles as described below.

4. Wash in a warm water/dishwashing detergent solution and rinse with clean, warm tap water. Then soak with a commercial disinfectant or one part vinegar (>=5% acetic acid concentration) to three parts warm water (131°F-149°F/55°C-65°C) solution for 60 minutes. Thoroughly rinse with warm tap water and air dry. **NOTE**– If using a commercial disinfectant, follow manufacturer's recommended instructions and dilution ratios carefully.

NOTE– Elbow is permanently attached to the lid and cannot be removed. Lid/Elbow assembly cannot be autoclaved.

Suction Unit (single-patient use)

1. With the power "Off," disconnect the unit from all external power sources.
2. Wipe the outside housing with a clean cloth dampened with commercial disinfectant/detergent (bacterial-germicidal).

CAUTION– Do not submerge suction unit in water; this will result in damage to vacuum pump.

6' Patient Tubing (single-patient use)

1. Disconnect the tubing from the unit.
2. Rinse thoroughly by running warm tap water through it.
3. Follow by soaking in a solution of 1 part vinegar (>=5% acetic acid concentration) to 3 parts water (131°F-149°F/55°C-65°C) for 60 minutes. Rinse with clean, warm water and air dry.
4. Keep the outer surface of the tubing clean by wiping with a clean, damp cloth.

Carrying Case (single-patient use)

1. Wipe the case using a clean cloth dampened with detergent and/or disinfectant.

NOTE– Disinfection information is based on AARC Clinical Practice Guideline Suctioning of the Patient in the Home.

Suction Unit (multi-patient use)

Device Cleaning and Disinfection When There is a Patient Change

DeVilbiss Healthcare recommends that at least the following procedures be carried out by the manufacturer or a qualified third party between uses by different patients.

NOTE– When the device is used as per instructions under normal conditions the interior of the device is protected from exposure to pathogens by the in-line filter on the collection container, hence no disinfection of internal components is necessary.

NOTE– If the device is used without an in-line filter then the interior of the device may be exposed to pathogens and the device cannot be disinfected.

NOTE– If the following described complete processing of the device by an appropriately trained individual is not possible, the device should not be used by another patient.

1. Use disinfectants safely. Always read the label and product information before use.
2. Always wear personal protective equipment when performing this procedure. Use suitable gloves and safety glasses. Cover exposed skin on arms to prevent accidental contact with bleach solution that has been applied to the device. Perform disinfection in a well-ventilated area.
3. Dispose of all accessories that are not suitable for reuse. This includes but may not be limited to the collection container, filter, tubing and carrying case.
4. With the power switch in the "Off" position, disconnect the device from all external power sources.
5. Visually inspect the device for any damage, missing parts, etc.
6. Clean the exterior of the device with a clean lint-free cloth. Heavy soil should be removed with a clean lint-free cloth dampened with water. A soft bristled brush dampened with water can be used to remove stubborn soil. Dry the device using a clean lint-free cloth if water was used to remove soil.
7. Use 5.25% chlorine bleach (Clorox Regular Liquid Bleach or equivalent). Mix one (1) part bleach with four (4) parts water in an appropriate clean container. This ratio produces a one (1) part bleach to five (5) total parts solution (1:5). The total volume (amount) of solution required is determined by the number of devices in need of disinfection.

NOTE– An alternate suitable commercial (bacterial-germicidal) disinfectant that meets the requirements listed in the NOTE below and is used as per the disinfectant manufacturer's recommended dilution rates and instructions may also be used.

CAUTION– Do not submerge the device in water as this will result in damage to the vacuum pump.

NOTE– Do not use any cleaners or disinfectants that contain ammonia, benzene and/or acetone to clean the device.

8. Apply the bleach solution in an even manner to the enclosure and power cord using a clean lint-free cloth. The cloth should be dampened only and not dripping of solution. Do not use a spray bottle to apply the solution. Do not saturate the device with the solution. Avoid over-saturating the enclosure seams so that no solution residue builds up in these areas.
9. Exposure time of the disinfectant solution should be 10 minutes minimum to 15 minutes maximum.
10. After the recommended exposure time, all exterior surfaces of the device are to be wiped with a clean lint-free cloth dampened with drinking quality water no warmer than room temperature. Dry the device with a dry, clean lint-free cloth. This is to remove residue that may stain or leave a film on the device, especially after repeated disinfections.

Disinfection

NOTE – The disinfection process can only be completed by the manufacturer or by an appropriately trained individual.

	Recommended disinfection interval	Number of disinfection cycles	Compatible disinfection method
Enclosure, power cord	Between patients	20	1:5 chlorine bleach (5.25%) and water solution, Mikrobac forte, Terralin Protect
Collection container, filter, tubing and carrying case	Do not clean, replace between patients	N/A	N/A

MAINTENANCE

Changing Filter & Replacing Container

1. Change the bacteria filter every 1 to 2 months OR immediately if overflow occurs.
2. Remove the bacteria filter by disconnecting it from the suction unit and lid assembly.
3. Replace it with a new DeVilbiss bacteria filter # 7305D-608 (12 pack).

NOTE– Verify clear side of filter marked <In> faces container.

NOTE– Use only the bacteria filter provided by DeVilbiss Healthcare or one of its Distributors. Substitution may lead to contamination of the unit and/or poor performance and will void warranty.

4. Replace Collection Container every 1 to 2 months after a minimum of 30 cleaning cycles as described in the cleaning instructions.

NOTE– Use only DeVilbiss 725 ml reusable containers. Substitution may lead to poor performance and will void warranty.

Service Life

- Service life of the suction unit is 5 years.
- Service life of the battery is 500 discharge cycles.
- Service life of the reusable container is 30 dishwasher cycles.

TROUBLESHOOTING

NOTE– Before you contact your equipment provider, follow Troubleshooting below.

DANGER

Electric shock hazard. Do not attempt to open or remove cabinet, there are no user-serviceable internal components. If service is required, return the suction unit to a qualified DeVilbiss provider or authorized service center. Opening or tampering with the unit will void the warranty.

NOTE– Your DeVilbiss Suction Unit contains no user-serviceable parts. If you believe your unit is not working properly, BEFORE YOU RETURN IT TO YOUR DEVILBISS HEALTHCARE PROVIDER WHERE YOU PURCHASED IT, PLEASE TAKE A FEW MOMENTS TO CHECK FOR THESE POSSIBLE CAUSES:

Problem	Action
Unit does not power on. (Green external power indicator should be illuminated when power is applied.)	<ol style="list-style-type: none">1. Check power sources and connections.2. Ensure wall outlet is live by plugging in a lamp.3. Check that battery is fully charged.
Pump runs, but no vacuum.	<ol style="list-style-type: none">1. Check that all tubing is connected properly.2. Check tubing connections for breaks, leaks, or occlusions.3. Ensure that float shut-off is not activated due to full container.4. Check for leaks or cracks in container assembly.
Low vacuum.	<ol style="list-style-type: none">1. Use vacuum adjustment knob to increase vacuum level.2. Check system for leaks.
Battery will not hold a charge.	<ol style="list-style-type: none">1. Confirm battery is connected (refer to Set-Up).2. Verify that charge light turns on.3. Check electrical connections during charging.4. Ensure wall outlet is live by plugging in a lamp. If problem continues, contact your equipment provider.

PROVIDER'S NOTES

No routine calibration or service is required provided the device is used in accordance with the manufacturer's directions. In case of a change of patient, the device must be reconditioned to protect the user. Reconditioning must only be carried out by the manufacturer or service provider. Between patients:

1. Visually inspect unit for any damage, missing parts etc.
2. Ensure that unit and accessories are clean.
3. Use an independent vacuum gauge to verify the unit provides the proper vacuum level as stated in Specifications.
4. Discard and replace collection container, filter, tubing, and carrying case between patients.
5. Wipe the surface using a clean cloth dampened with disinfectant.

SPECIFICATIONS/CLASSIFICATIONS

Size	7.25"H x 7.25"W x 6.75"D (18.4cm x 18.4cm x 17.1cm)
Weight	3.375 lbs. (1.53 kg)
Electrical Requirements	100-240 VAC 50/60 Hz 0.75 A max; 12 VDC, 33 W max
Vacuum Range	50-500 mmHG +/- 10%
Air Flow @ pump inlet	27 LPM (free flow) typical (may be less when running from internal battery)
Collection Container Capacity	725 ml (cc)

Environmental Conditions

Operating Temperature Range	32°F (0°C) - 104°F (40°C)
Operating Relative Humidity	0-95%
Operating Atmospheric Pressure	10.2 Psi (70 kPa) - 15.4 Psi (106 kPa)
Storage & Transport Temperature Range	-40°F (-40°C) - 158°F (70°C)
Storage & Transport Relative Humidity	0-95%
Storage & Transport Atmospheric Pressure	7.3 Psi (50 kPa) - 15.4 Psi (106 kPa)

Warranty

7310 Series	Two-years limited, excluding internal battery and collection container
Internal Battery	.6 months

Approvals

IEC 60601-1:2005+AMD1:2012; EN 60601-1-2: 2015; ISO 10079-1:2015; IEC 60601-1-2:2014; IEC 60529:1989+AMD1:1999 +AMD2:2013; IEC 60601-1-6:2010; IEC 60601-1-11:2015; IP22
CONFORMS TO AAMI STD. ES 60601-1, AAMI STD. HA 60601-1-11, IEC STD. 60601-1-6, ISO STD. 10079-1
CERTIFIED TO CSA STD. C22.2 NO. 60601-1, NO. 60601-1-11, IEC STD. 60601-1-6, ISO STD. 10079-1
7310 Series complies with IEC 60601-1 3rd Edition

Equipment Classifications

With respect to protection from electric shock Class II and internally powered
Degree of protection against electric shock Type BF Applied Parts
Degree of protection against ingress of liquids IP22 and ordinary power supply
Mode of Operation Intermittent Operation: 30 minutes on, 30 minutes off
Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

ISO Classification

7310 Series - Electrically powered medical suction equipment for field and transport use according to EN ISO 10079-1:2015
High Flow/High Vacuum

TWO-YEAR LIMITED WARRANTY

The compressor portion of the DeVilbiss Vacu-Aide Compact Suction Unit 7310 Series (excluding internal rechargeable battery and collection container) is warranted to be free from defective workmanship and materials for a period of two years from date of purchase. Internal rechargeable batteries are warranted for 6 months. Any defective part(s) will be repaired or replaced at DeVilbiss Healthcare's option if the unit has not been tampered with or used improperly during that period. Make certain that any malfunction is not due to inadequate cleaning or failure to follow the instructions. If repair is necessary, contact your DeVilbiss Provider or DeVilbiss for instructions.

NOTE– Be sure to retain a dated proof of purchase document to verify unit is within 2-year warranty period.

NOTE– This warranty does not cover providing a loaner unit, compensating for costs incurred in rental while said unit is under repair, or costs for labor incurred in repairing or replacing defective part(s).

THERE IS NO OTHER EXPRESS WARRANTY. IMPLIED WARRANTIES, INCLUDING THOSE OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED TO THE DURATION OF THE EXPRESS LIMITED WARRANTY AND TO THE EXTENT PERMITTED BY LAW. ANY AND ALL IMPLIED WARRANTIES ARE EXCLUDED. THIS IS THE EXCLUSIVE REMEDY AND LIABILITY FOR CONSEQUENTIAL AND INCIDENTAL DAMAGES UNDER ANY AND ALL WARRANTIES ARE EXCLUDED TO THE EXTENT EXCLUSION IS PERMITTED BY LAW. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, OR THE LIMITATION OR EXCLUSION OF CONSEQUENTIAL OR INCIDENTAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state.

Manufacturer's Note

Thank you for choosing a DeVilbiss Vacu-Aide Compact Suction Unit. We want you to be a satisfied customer. If you have any questions or comments, please send them to our address on the back cover.

For Service Call Your Authorized DeVilbiss Healthcare Provider:

Phone	Purchase Date	Serial #
-------	---------------	----------

ELECTROMAGNETIC COMPATIBILITY INFORMATION

WARNING

MR Unsafe

Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the Vacu-Aide or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.

Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the Vacu-Aide. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.

This device is suitable for use in home and healthcare environments except for near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of Electromagnetic DISTURBANCES is high.

Cables and maximum lengths of cables

DC power cord (cigarette lighter adapter) #7304D-619 maximum length = 2.2 meter (7.2 ft)

WARNING

Use of accessories and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.



WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Vacu-Aide, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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DeVilbiss® and Vacu-Aide® are registered trademarks of DeVilbiss Healthcare.
© 2019 DeVilbiss Healthcare LLC. 12.19 All Rights Reserved.
Drive DeVilbiss Healthcare Ltd.
Drive DeVilbiss Australia Pty Limited
DeVilbiss Healthcare GmbH
DeVilbiss Healthcare LLC • 100 DeVilbiss Drive • Somerset, PA 15501 • USA