

*drive*

 DeVilbiss<sup>®</sup>  
HEALTHCARE

# POWERCHAIRS

PRACTICAL. ELECTRIC. UNLIMITED MOBILITY.

# SIMPLY. INDEPENDENT. MOBILE.

Limited mobility does not mean a limited lifestyle. Our power wheelchairs make sure of that. Drive DeVilbiss Healthcare lightweight folding power wheelchairs are ideal for supporting people with limited mobility and those who struggle to walk long distances.

## SIMPLY. INDEPENDENT. MOBILE.

For those who like to keep on the move!

For everyday life or even when travelling, an electric power wheelchair is the perfect companion. Once folded, these power wheelchairs are on board in no time at all, whether you are travelling by car, train, ship or even on an aeroplane\*, the Instafold and Airfold power wheelchairs give you total independence and the freedom to explore the world for what it has to offer.

\*Always clarify the conditions of carriage with the airline before travelling.



**Intelligent  
Control**



**Manual  
folding**



**Lithium-ion  
Battery**

# INSTAFOLD



## Our entry-level aluminum model

Folded in an instant and ready to go! The **Instafold** power wheelchair, made from aluminum, not only offers stability and durability, this economical powerchair solution impresses with its modern stylish design.

The Instafold is compact, **weighs only 27 kg including the 20Ah lithium-ion battery** and provides maximum comfort for users up to 120 kg body weight. With an intelligent control feature the powerchair is user friendly, easy to operate and promotes unlimited independence.

# AIRFOLD



## Our Premium lightweight model

The **Airfold** powerchair is a new design of innovation, constructed from corrosion resistant carbon fibre material with a trendy elegant design.

This premium compact folding powerchair is ultra-lightweight and **weighs only 18.5 kg, including its 20Ah lithium-ion battery** and provides a maximum user weight of 115 kg. Supplied in a trolley style storage bag with integrated handles, the Airfold can be easily transported and stored when not in use.

## INSTAFOLD

Our entry-level aluminum model



Item code

IFPC17SIL-FR

## AIRFOLD

Our Premium lightweight model



Item code

ACFPC17BLK-EU

## TECHNICAL DATA

Speed	6 km/h
Battery Capacity	1 x 20 Ah lithium Ion
Turning Radius	196 cm
Max. Gradient	10 %
Ground Clearance	7,5 cm
Measurement (LxBxH)	105 x 62 x 99 cm
Folding Measurements (LxBxH)	79 x 62 x 40 cm
Seat Width/Depth	44 cm / 45 cm
Product Weight w/o Battery	24,3 kg
Max. User Weight	120 kg

Speed	6 km/h
Battery Capacity	1 x 20 Ah lithium Ion
Turning Radius	120 cm
Max. Gradient	10 %
Ground Clearance	6,5 cm
Measurement (LxBxH)	93 x 57 x 93 cm
Folding Measurements (LxBxH)	79 x 57 x 32 cm
Seat Width/Depth	44 cm / 42 cm
Product Weight w/o Battery	15,8 kg
Max. User Weight	115 kg

# DRIVE DEVILBISS, A POWER BRAND YOU CAN TRUST.

We are a Global leading manufacturer and supplier of durable medical equipment with a long-standing company history.

With a successful range of innovative mobility products and electric scooters, we have applied that knowledge and know-how and developed a new range of electric powerchairs. **Why not see for yourself!**



Drive Medical GmbH & Co. KG • Leutkircher Straße 44 • 88316 Isny/Allgäu  
DeVilbiss Healthcare GmbH • Kamenzer Straße 3 • 68309 Mannheim  
T +49 (0) 7562 9724 0 • info@drivemedical.de • [www.drivedevilbiss.de](http://www.drivedevilbiss.de)