

# Product Information – Technical Change

Dear customers and business partners,

Due to constant changes within the market place, and to ensure our products meet market requirements, we would like to inform you about a technical change regarding our Lightweight Rollator Gigo.

As of today, the 'Lightweight Rollator Gigo 2G' will be a direct replacement for the Rollator Gigo. This product development includes the following technical changes:

- *Removable and folding backrest, for more flexibility and ease of use when folding the rollator*
- *Adding a Curb Climber as standard*

Please see item number(s) details below:

## Replacement product:

Lightweight Rollator Gigo 2G

Item number: 735100200

HMV-number: 10.50.041256

## Discontinued product:

Lightweighthrollator Gigo

Item number: 735100100

No longer available from: immediately



For more information, please visit our website: [www.drivedevilbiss.de](http://www.drivedevilbiss.de)

Should you have any further questions about this notice or other products, please contact your local Drive DeVilbiss representative.

Kind regards

Your Drive DeVilbiss Team