


# Overbed table SuppliTo

For support in daily life



Special angle fixation



The SuppliTo overbed table is intended for use in the home environment and is designed to provide support for homecare. SuppliTo supports the user and patient at the same time and ensures a very high level of patient comfort thanks to various adjustment options. The table top can be adjusted in inclination and offers maximum safety thanks to special fixations.

- Continuous height adjustment
- Virtually tool-free disassembly
- 4 smooth-running rollers, 2 with brakes
- Flexible use in the patient room
- Individual inclination adjustment of the table top with inclination fixation

Product	white-grey	wood optic
Order number	860000001	860000003

Technical Data	
Overall height	71 - 116 cm
Length tabletop	60 cm
Width tabletop	40 cm
tilt adjustment	0 - 90 °
Total weight	8 kg
Max. Load	10 kg

# Overbed table SuppliTo twinTop

For support in daily life



The SuppliTo twinTop overbed table consists of a tilt-adjustable and a fixed tabletop. A special inclination fixation ensures maximum safety. It supports both the user and the patient, e.g., during homecare activities or food intake in bed. The height of the SuppliTo twinTop can be individually adjusted. Even at higher table settings, the overbed table remains stable.

- Easy to fold up
- Tool-free assembly and disassembly
- 4 smooth-running rollers, 2 with brakes
- Additional table top
- Flexible use in the patient room

Product	white-grey	wood optic
Order number	860000002	860000004

Technical Data	
Overall height	71 - 93 cm
Length tabletop	60 und 20 cm
Width tabletop	40 cm
tilt adjustment	0 - 50° in 7 Schritten
Total weight	9 kg
Max. Load	10 kg